

# Sklar Disinfectant

**3-Minute Ready-To-Use Disinfectant, Decontaminant & Cleaner  
For Medical, Dental, Ophthalmic and Surgical Equipment & Surfaces**

**Bactericidal – Virucidal – Fungicidal – Tuberculocidal**

**DESCRIPTION:** Sklar Disinfectant is a ready to use, highly effective cleaner and disinfectant for use in intensive care units and high-risk patient care areas, surgical suites, infant care units, plasma phoresis facilities and other areas where strict environmental control of cross contamination is critical.

- Sklar Disinfectant's non-caustic formula is used to thoroughly clean and disinfect dental, ophthalmic, medical and surgical equipment surfaces such as operating room lights and tables, IV and chest pumps, EMS equipment, examination room fixtures, infant incubators, respiratory therapy and anesthesia equipment, laboratory equipment and surfaces as well as other contaminated items.

- Sklar Disinfectant may be applied directly to surfaces or used as an immersion solution. Complete instructions for its variety of uses are listed on its label.



**MICROBIOLOGY:** Sklar Disinfectant completely kills infectious disease causing microorganisms within 3 minutes, including:

**Gram Positive and Negative Bacteria**

Mycobacterium tuberculosis (TB)  
Staphylococcus aureus  
Pseudomonas aeruginosa  
Salmonella chloreraesuis

**Viruses**

Hepatitis B Virus (HBV)  
Hepatitis C Virus (HCV)  
Influenza A/2 Hong Kong Virus  
Rhinovirus  
Human Immunodeficiency Virus (HIV-!)  
Herpes Simples Virus 1 and 2  
Polio Virus

**Pathogenic Fungi**

Trichophyton mentagrophytes

- **Broad Spectrum Ready-To-Use Disinfectant/Cleaner**
- **Kills TB/Hepatitis B & C In 3 Minutes!**



**SAFETY:** Sklar Disinfectant is safe for use on inanimate surfaces including plastics, stainless steel, chrome, glass, Plexiglas® and baked-on painted surfaces. Sklar Disinfectant leaves no residue.

Sklar Disinfectant meets OSHA blood-borne pathogen requirements for disinfection in healthcare facilities and is formulated following strict FDA Good Manufacturing Practices.

**APPLICATIONS:** When used according to label instructions, Sklar Disinfectant is highly effective:

- Surface Disinfectant and Cleaner
- Instrument Decontaminant Solution
- Ultrasonic Disinfection Cleaning Solution

<u>Ordering Information</u>	<u>Product No.</u>	<u>Packaged</u>
24 oz Bottles	#1643	Each
One Gallon Bottle	#1653	Each

Special Case (of 4 gallons) Pricing Available

## Dental Corporation of America

889 South Matlack Street  
West Chester, PA 19382

Phone: 800-638-6684 • 610-344-7488

Fax: 610-431-6500

ortho@dentalcorp.com

[www.dentalcorp.com](http://www.dentalcorp.com)